


## Cardiac Troponin I Rapid Test Cassette (Whole Blood /Serum/Plasma)


### Whole Blood specimen:

- transfer 3 drops (approximately 75 µl) of fingerstick whole blood specimen to the specimen well (S) of test cassette, then add 1 drop of buffer (approximately 40 µl) and start the timer.



### For Serum or Plasma specimen:

- transfer 3 drop of serum or plasma (approximately 75 µl) to the specimen well (S) of test cassette, then add 1 drop of buffer (approximately 40 µl) and start the timer.



\*Note: Wait for the colored line(s) to appear. **Read results at 10 minutes.** Do not interpret the result after 20 minutes.


### INTERPRETATION OF RESULTS

**POSITIVE:** \* **Two lines appear.** One colored line should be in the control line region (C) and another apparent colored line should be in the test line region (T).

**\*NOTE:** The intensity of the color in the test line region (T) will vary depending on the concentration of cardiac Troponin I (cTnI) present in the specimen. Therefore, any shade of color in the test line region (T) should be considered positive.

**NEGATIVE:** **One colored line appears in the control line region (C).** No line appears in the test line region (T).

**INVALID:** **Control line fails to appear.** Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit immediately and contact your local distributor.



### PRECAUTIONS

- This package insert must be read completely before performing the test. Failure to follow directions in insert may yield inaccurate test results.
- For professional in vitro diagnostic use only.
- Do not use after expiration date.
- Wear protective clothing such as laboratory coats, disposable gloves, and eye protection when specimens are assayed.
- Please ensure that an appropriate quantity of sample volumes are used for testing. Too much or too little sample size may lead to deviation of results.
- The used test should be discarded according to local regulations.

### STORAGE AND STABILITY

Store as packaged in the sealed pouch at room temperature or refrigerated (2-30°C). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use.

**DO NOT FREEZE.** Do not use beyond the expiration date.